## **ACTREC**

## TEST DIRECTORY REPORT LOCATION DETAILS

Department name : CANCER CYTOGENETICS DEPARTMENT

Location : Room No. 6, Ground Floor, CCE Building.

## **FUNCTIONAL DETAILS**

Test Name : Constitutional karyotyping

Test code : T311
Test status : ACTIVE
Modality : CK

Quantity required : Peripheral Blood: 5 ml

Collection Instructions : Refer Primary sample collection manual

Nature of specimen : Peripheral Blood

Method used : Conventional Karyotyping

Type of container : Sterile sodium heparin Green top vacutainer
Transport instructions : Refer Primary sample collection manual

Time for additional exam : NA

on stored sample, if applicable

Resource person for report status : Officer-in-charge Resource person for test query : Officer-in-charge

Patient instruction : Refer Primary sample collection manual Dept. acceptance time : Monday to Friday-9:30 a.m. till 5:00 p.m.

Saturday- 9:30 a.m. till 12:00 p.m.

Turn around time : 17 working days

Sample storage after reporting : Fixed pellet stored in minus 20 deep freezer

Telephone : 022-27405000

Test schedule : Daily

Extension : 5758, 5759

Use : Constitutional abnormalities