ACTREC **TEST DIRECTORY REPORT** LOCATION DETAILS

Department name : CANCER CYTOGENETICS DEPARTMENT : Room No. 6, Ground Floor, CCE Building.

Location

FUNCTIONAL DETAILS

Test News	· A suto Louisonaio komuntumine
Test Name	: Acute Leukemia karyotyping
Test code	: T315
Test status	: ACTIVE
Modality	: CK
Quantity required	: Bone Marrow: 1-4 ml, Peripheral Blood: 4-5 ml
Collection Instructions	: Refer Primary sample collection manual
Nature of specimen	: Bone marrow, Peripheral blood (Presence of more than
	80% blasts in peripheral blood)
Method used	: Conventional Karyotyping
Type of container	: Sterile sodium heparin Green top vacutainer
Transport instructions	: Refer Primary sample collection manual
Time for additional exam	: NA
on stored sample, if applicable	
Resource person for report status	: Officer-in-charge
Resource person for test query	: Officer-in-charge
Patient instruction	: Refer Primary sample collection manual
Dept. acceptance time	: Monday to Friday-9:30 a.m. till 5:00 p.m.
	Saturday- 9:30 a.m. till 12:00 p.m.
Turn around time	: 17 working days after hematopathology report
Sample storage after reporting	: Fixed pellet stored in minus 20 deep freezer
Telephone	: 022-27405000
Test schedule	: Daily
Extension	: 5758, 5759
Use	: Acute leukemia